
900-38

900-39

(EN) User manual bathroom scale

(EN) 1. Specifications

- 1-1. Touch keys: SET / (up) / (down)
- 1-2. Batteries: 3 × AAA
- 1-3. Maximum load: 180 kg / 397 lb / 28 st: 5 lb
- 1-4. Graduation: 0.1 kg / 0.2 lb / 1 lb
- 1-5. Overload indication: Err
- 1-6. Low battery indication: Lo
- 1-7. Zero indication: - - -
- 1-8. Height range: 100 ~ 240 cm / 3ft 3inch ~ 7ft 10inch
- 1-9. () Age range: 10-80 years old
- 1-10. () Number of users: 10 persons (P0-P9)
- 1-11. () Measuring body fat range: 5 % ~ 50 %
- 1-12. () Measuring water range: 10 % ~ 80 %
- 1-13. () Measuring muscle range: 10 % ~ 80 %
- 1-14. () Measuring bone range: 0.5 kg ~ 8.0 kg
- 1-15. () Recommended daily intake of calories
- 1-16. () BMI (Body Mass Index)
- 1-17. Automatic user identification range: 2.0 kg
- 1-18. Normal mode and Sport mode ()

(EN) 2. Important!

- 2-1. Ensure the scale is on a flat, even and firm surface. Avoid carpet or soft surfaces.
- 2-2. The surface of the scale may become slippery if it gets wet. Please make sure that both the glass and your feet are dry before using the scale.
- 2-3. This scale is not water proof. Do not expose it to moisture.
- 2-4. Clean the scale with a slightly damp cloth. Do not use any chemical/abrasive cleaning agents to clean the scale.
- 2-5. Remove the batteries if the scale is supposed to be used for a longer time.
- 2-6. The appliance is intended solely for domestic use within the home. It is not intended for professional use.
- 2-7. Children should be supervised to ensure that they do not play with the appliance.
- 2-8. Before first time use: Carefully unpack the appliance and remove all the packaging material. Keep the material (plastic bag and cardboard) out of reach of children. Remove the batteries from the protective packaging and insert them with the correct polarity.

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
(EN) 3. Unit

3-1. Press the UNIT button on the back of the scale to select measuring unit (kg/lb/st). The default unit is kg.

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(EN) 4. User parameters

You must enter user parameters in order to measure your body fat percentage and other physical data, by following the steps below.

Memory position (P)	P0-P9
Body height	100~240cm 3ft 3inch ~ 7ft 10inch
Age	10-80 years
Gender	Male (♂), female (♀)
Sport mode	Sport: Male (♂), female (♀) For people who exercise regularly

4-1. When the scale's display is active, press SET to enter the setting mode.

4-2. Press or to choose a user number. Press SET to confirm.

4-3. Press or to choose Male, Female, Sport-Male or Sport-Female (use one of the two Sport modes if you exercise regularly). Press SET to confirm.

4-4. Press or to enter body height. Press SET to confirm.

4-5. Press or to change the age. Press SET to confirm.

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
(EN) 5. Weighing and measuring body data (first time use)

5-1. When the user parameters are set, the display will show "0.0".

5-2. Stand barefoot on the scale and make sure that your feet have good contact with the electrode plates in metal. Your weight will show up in a second.

5-3. Stay on the scale while it is measuring your body data. Let the progress bar finish ("OOOO").

5-4. Your body data will be displayed one by one: Body fat ratio, body water ratio, BMI, muscle mass, bone mass, calories and weight. All body data will be displayed twice.

5-5. Except user number and gender, your height and age will also be shown for you so that you can confirm them. The scale will eventually turn itself off.

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

(EN) 6. Automatic user identification (range 2.0 kg)

The next time you use the scale, it will automatically identify you based on your weight if the weight is within a range of 2.0 kg and not in conflict with another user's weight. It means that you can just stand on the scale without having to select your user number or re-enter the user parameters.

6-1. Stand barefoot on the scale and make sure that your feet have good contact with the electrode plates in metal. Your weight will show up in a second.

6-2. Stay on the scale while it is measuring your body data. Let the progress bar finish ("OOOO").

6-3. Your body data will be displayed one by one: Body fat ratio, body water ratio, BMI, muscle mass, bone mass, calories and weight. All body data will be displayed twice.

6-4. If the scale can't identify your user number, e.g. two user's weight are too close to each other, you should press  or  to select your user number.

6-5. If the scale is only showing your weight and not showing any other body data, then make sure that you are barefoot and that your feet have good contact with the electrode plates in metal.

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
(EN) 7. Weighing only

7-1. If you only want to measure your weight, step onto the scale and wait. Your weight will appear on the display in a second. Step off the scale. It will automatically turn off in a few seconds.

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(EN) 8. Evaluation of the results

The body data tables below are guidelines only.

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(EN) According to human body basic metabolism, the body needs a certain amount of calories every day. The value that is shown on the scale display and in the table below is a recommended daily intake of calories for people who are physically active, e.g. using weights.

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